

# Adaptive Stress Reduction Module

What it is:

- A chest worn module to track physiological state paired with a breathing software designed to give personalized breathing routines to reduce stress

What it does:

1. Measures users stress via chest worn sensor array
2. Use ML/LLM to prescribe personalized breathing exercise
3. Continuously adjust meditation guidance to optimize stress reduction based on biometric feedback

Key Features


1. Real time biometric monitoring (ECG, breathing movement, temperature)
2. AI driven meditation practices with real time adjustments personalized to the user
3. Session feedback and stress reduction metrics

Workflow:



1. Session Initiation → Baseline Measurement → AI Prescription
2. Active Meditation → Continuous Monitoring → AI driven Adjustment to optimize session
3. Session Completion → Analysis → Results presentation

Data Flow

User Starts Session



Raw Sensor Data



MTTQ Publish  
Base line Data

Processing System  
Laptop/Device/App

User Input

- Stress Level
- Time available

Data Processing

Biometric Data → Stress Score/Category  
Meaningful Context

LLM Analysis  
(Chat GPT/Claude/Gemini)

Input/Prompt

- Stress Assessment
- Category
- Primary Indicators
- User Input

LLM

Response/Meditation Prescription

- Pattern
- Inhale:
- Hold:
- Exhale

Active Meditation

Wearable Sensor Module

Continuous Reading (100ms?)

- ECG Waveform
- Chest Movement
- Skin temp

Real time Processing

- R-R intervals
- Breathing Cycles
- Pattern Adherence

MTTQ  
Publisher → Subscriber

Device/App

- Initiate breathing pattern and associated interface

Real time Adjustments  
(Every Xsec or Ymin)

- Store the data (JSON)
- Process the data Frame
- Prompt LLM
  - Should any adjustments be made
- LLM Response
  - Adjustment Command
- Silently update interface

Repeat till  
Duration is  
Finished

LLM

Completion/Exit

### Final Metrics (Wearable)

#### Continuous Reading

- ECG Waveform
- Chest Movement
- Skin temp

- R-R intervals
- Breathing Cycles
- Pattern Adherence

MTQ  
Publisher → Subscriber

### Device/App

- Store Data (JSON)
- Session analysis

Note: Data Storage could be

- SQLite
- JSON

Would Store


- User Baseline
- Session history
- Effectiveness metrics

### Visualization Dashboard

#### Stress Reduction Achieved


#### Heart Rate

Average: 74 bpm



#### HRV

Average: 52 66 •



#### AI Insights

"Excellent! .....

### Store the Session

- Store key metrics of session in SQLite
  - User Baseline
  - Session Info
  - Effectiveness Metric

